

Clip Finger Pulse Oximeter

夾指式脈搏血氧儀

Model : HM-DV

Features:

- Small and light, one key measurement, simple operation, use more worry
- 130 degree viewing angle, TFT LCD color screen, screen brightness 4 degrees adjustable
- Silicone finger mold, soft and comfortable, does not press the finger tip, blood flow is more smooth
- 4-direction display, clearer reading, press the switch button to adjust
- It is suitable for family, oxygen bar, exercise health care (before and after exercise, not during exercise)
- 小巧輕便 一鍵測量, 操作簡單, 使用更省心
- 130度可視角, TFT液晶彩屏, 螢幕亮度4度可調
- 矽膠指模, 柔軟舒適, 不壓迫手指末梢, 血流更流暢
- 4方向顯示, 讀數更清晰, 輕按開關按鍵調節
- 適用於家庭、氧吧、運動保健(在運動前後使用, 不建議在運動過程中使用)等範圍



Specifications

Product Model 產品型號	HM-DV
Product Function 產品功能	blood oxygen, heart rate and pulse test 血氧、心率脈搏測試
Product Material 產品材質	ABS, TFT LCD display ABS, TFT液晶屏顯示
Measuring Range 血氧測量範圍(SPO ₂)	100% ≤ normal ≤ 94% ≤ insufficient oxygen supply 70~100% (±2%)
脈搏測量範圍(PR)	60bpm ≤ normal ≤ 90bpm 30~240bpm (±2bpm)
血流灌注指數(PI)	The greater the pulsating blood flow, the greater the PI value 脈動的血流越大, PI值就越大
Power Supply 電源方式	AAA batteries *2 AAA電池*2
Product Size 產品尺寸	5.7*3.3*3.2cm
Product Weight 產品重量	27g
Packing Size 包裝尺寸	9*6*3.8cm
Package Weight 包裝重量	70g

Instructions

- ◆ 一鍵開機: 一鍵開機功能, 操作簡單
- ◆ 正確測量是將指甲蓋這面向上放置於血氧儀的測試位置, 打開開關血壓儀將自動測試并顯示結果
- ◆ 開機後, 長按開關機3-5秒可進入設置菜單, 找到你需要設置的菜單長按3秒左右即可自己設置

注意事項:

- ◆ 當你的手指進入血氧儀時, 你的指甲表面必須向上(每次測試前請用酒精清洗橡膠, 測試前後請用酒精清洗被測手指)
- ◆ 此產品不適合作為病人的持續監測使用
- ◆ 當電量低時應及時換上電池, 避免影響正常使用